

Standard Operating Procedure for Hazardous Chemicals

Principal Investigators: Chung-Jui Tsai and Scott A. Harding

Building and rooms: Davison Life Sciences Building, Lab B310

Chemical(s)	Isopropanol, chloroform, ethanol, liquid nitrogen (LN), SDS, ethylenediamine tetraacetic acid (EDTA), isoamyl alcohol (IAA), Itris(hydroxymethyl)aminomethane (Tris-base), sodium hydroxide (NaOH), potassium acetate.
Process	TA cloning
Specific Hazards <i>referred to MSDSs for more detailed information</i>	<u>Chloroform</u> is a probable human carcinogen, reproductive hazard and skin/eye irritant. <u>IAA</u> is flammable, and may be harmful upon skin/eye contact, ingestion or inhalation. <u>LN</u> is a frostbite hazard. <u>NaOH</u> is a corrosive base, may cause burns to skin and eyes upon contact and to mucous membranes if inhaled or ingested.
Personal protective equipment	Must wear 3-5 mil nitrile gloves. Chemical safety goggles and lab coat should be worn when splash potential exist
Engineering/ventilation controls	All operations involving <u>chloroform</u> , and <u>IAA</u> must be done in a chemical fume hood.
Special handling procedures and storage requirements	Store <u>chloroform</u> , <u>IAA</u> , and <u>alcohol</u> in the flammable cabinets under the hood. Store all the reagents on bench top, 4°C or -20°C accordingly.
Spill and accident procedures <i>for hazardous chemicals only</i>	<p><u>Skin exposure</u>: Rinse affected skin with plenty of water while removing contaminated clothing/shoes. Rinse for > 15 minutes.</p> <p><u>Eye exposure</u>: Wash eyes for > 15 minutes.</p> <p>For both cases, seek medical attention immediately.</p> <p><u>Small</u> (< 2L): Absorb with vermiculite or spill pads and transfer absorbed material to a closed container. Label and date as hazardous waste for disposal. Notify PIs.</p> <p><u>Large</u> (> 2L): Evacuate the room, notify PIs and call 2-5801 to request emergency spill assistance from the Environmental Safety Division.</p>
Waste disposal	Chloroform, IAA waste must be collected and labeled as hazardous waste according to the SOP for Hazardous Waste Disposal.
Special approval	No special authorization needed after SOP training & reading MSDSs.
Prepared by	Name/date: Kate Tay, 9/8/2010
Reviewed by	Name/date: C-J Tsai, 9/9/2010

TA Cloning

TA ligation (Invitrogen 461370, Dual Promoter TA Cloning Kit) or home-made pCX T-vectors
Ligation should be done **right after** PCR reaction. The A-overhangs on the PCR products are unstable overtime.

Half-reaction	needed	final	
10x Ligation buffer ¹	1 μ l	1 x	} Mix, short spin 14°C, O/N Stored at -20°C (can be re-used if needed, label clearly w/ vector + insert and date) OK to use NEB ligase as well
10mM ATP	0.5 μ l	0.5 mM	
pCRII or pCX vector	0.5 μ l	~10 ng	
Insert (fresh PCR product) ²	1 μ l	see below	
Ligase	0.5 μ l	2 units	
ddH ₂ O	<u>6.5μl</u> <u>10.0μl</u>		

¹ The 10X ligase buffer precipitates after thawing on ice. Vortex vigorously and spin down before use. NEVER thaw the ligase buffer at room temp.

² Shoot for a molar ratio of insert : vector > 3:1 ($[\text{insert length bp}]/[\text{vector length bp}] \times 10 \text{ ng} \times 3$).
For good PCR products, 1 μ l is usually sufficient. Too much DNA will reduce transformation efficiency.

Prepared LB+Kan plate

	Final	1L
LB (low salt) Powder	--	20 g
Agar	15%	15 g
[100mg/ml] Kan	100 mg/L	1 ml
[100mg/ml] Amp	100 mg/L	1 ml

} Autoclave & cool to touch
+ Antibiotic & pour

} 25 ml per plate;
40 plates per liter

Transformation

1. Turn on water bath to 42°C (or warm a beaker of water to 42°C by microwave)

2. Ligation mix 2 μ l } Tap only once to mix GENTLY!
Top10 C.C. (purple)* 15 μ l } Do Not pipette up & down!
On ice, 30 min

3. Heat Shock: 42°C, 30sec (Do Not mix or shake!) & chill on ice

4. + 500 μ l SOC/LB (in laminar flow)

5. Shake, 37°C, 1 hr, 225 rpm (Warm up LB + Kan plate at 37°C, 10min)

6. Take 200 μ l of the transformant mix

→ pCRII spread 200 μ l is more than enough

→ pCXSN series need to spread more, ~300 μ l

and + 20 μ l [40 mg/ml] X-gal (light sensitive)
+ 20 μ l [10 mM] IPTG

} **Double check!**
TOP 10 (general use) requires only X-gal
TOP 10F' (methylated/toxic genes) requires
both X-gal & IPTG
pCX series does not require any

5. Spread all on plate & dry

6. Incubate @ 37°C, O/N

- * Competent cells are expensive. They should be kept cold at all time and thaw only on ice.
- * Each vial has 50 μ l and can be used for up to 3 reactions (at least 15 μ l each).
- * For multiple transformations, add 2 μ l ligation reaction to a fresh tube on ice. Transfer 15 μ l of cells to each tube, then add the 3rd ligation to the remaining cells in the original CC tube.
- * If you have only 1 transformation to perform, use the whole vial or wait until you have more.
- * Re-frozen competent cells have drastically reduced transformation efficiency.

Colony PCR


For blue/white selection, incubate the plate at 4°C for > 1 hr

Colony Lysis buffer	stock	20 colonies
1% Triton X-100	100%	<u>10ml</u> 100 µl
20 mM Tris, pH8.5	500 mM	400 µl
2 mM EDTA	100 mM	200 µl
s/w		9.3 ml

1. Aliquot 50 µl of lysis buffer into each tube.
2. Number the (white) colonies. Pick a part of the colony with a pipet tip and mix with lysis buffer by pipetting up and down several times. Re-incubate the plate at 37°C for a few hrs.
3. 95°C, 10 min
4. Centrifuge: 2 min, RT, max. speed
5. Transfer supernatant to new tube (optional)
6. Take 1 µl lysate for colony PCR. Stored the rest at -20°C

7. Prepare **master mix** [# of reactions plus ~10% (or one reaction) extra] for PCR:

Lysate	1	µl
10x REDTaq buffer	1	µl
10mM dNTP	0.25	µl
5µM forward primer*	0.5	µl
5µM reverse primer *	0.5	µl
0.5 U REDTaq	0.5	µl
1.5mM MgCl ₂	3	µl
s/w	<u>3.25</u>	<u>µl</u>
	<u>10</u>	<u>µl</u>



* Use both gene-specific F + vector F AND gene-specific F + vector R to determine the insert orientation in the vector (expect one positive and one negative reaction).

8. Run 5µl on a gel.

For PCR-positive colonies:

Pick colonies (choose 4)

1. Prepare 5 ml LB liq + 5 µl [100 mg/ml] Kan in culture tubes (final [Kan] = 100 mg/L)
2. Wipe the pipet shaft with 70% EtOH. Pick the desired colony with a pipet tip and eject the tip into the LB media
3. Shake, 37°C, O/N, 225 rpm (8-16 hrs, the shorter the better, no more than 16 hrs!)
4. Use immediately or store at 4°C overnight for miniprep (or up to 2 weeks max for subculturing).

Glycerol Stock preparation:

1. Take 0.5 ml cultures for glycerol stock in a 2-ml cryo-tube + 0.5 ml LB (no antibiotic) w/ 80% glycerol (1:1 ratio)
2. Snap-freeze in LN before storing at -80°C.

Miniprep by alkaline lysis

GTE buffer (50 mM glucose, 25 mM Tris-HCl pH8, 10 mM EDTA pH8) Sol. 1 → on ice from 4°C
0.2M NaOH/1%SDS (make **fresh weekly** w/ 10M NaOH and 20% SDS) Sol. 2 → RT
5 M potassium acetate (pH4.8) Sol. 3 → on ice from 4°C
RNase A (DNase-free), 10 mg/ml
Chloroform:IAA = 24:1
Isopropanol

1. Transfer 1.5 ml cultures to a new 1.7 ml tube, centrifuge at RT, max. speed, 45 sec.
Discard supernatant directly. Repeat for up to 2 more times (in the same tube). Unused bacterial waste needs to be autoclaved.
2. Resuspend the pellet in 200µl of **cold** GTE (1) by pipetting up/down.
Disturb the pellet to ensure **complete suspension!** (Turns turbid).
3. Denature with 400µl of 0.2N NaOH/1% SDS (2), **gently** invert 5 times (become clear). NO incubation.
4. Neutralize with 300µl of **cold** 5M KOAc (3), **gently** invert (white protein precipitate form). NO incubation.
5. Centrifuge, RT, max, 10 min. Transfer supernatant (900µl clear, only DNA/RNA) to new 1.5ml.
6. Add 1.8 µl RNase A (invert gently once) to 20 µg/ml. Incubate at 37°C, 20-40 min.
7. Extract with 400 µl chloroform:IAA, Vortex 5 sec. Centrifuge: RT, max, 1 min.
Transfer top phase to new 1.5 ml (4x175 µl, tilt the tube to form a ball).
DO NOT touch the interphase or the bottom layer! Leave behind some top phase (quality is more important – you will get more than you need)
8. (Repeat once more #7). Chloroform waste goes into a specifically labeled waste jar.
9. Precipitate DNA with 1 Vol (700 µl) isopropanol, invert 5 times, and **immediately** centrifuge, RT, max, 10 min. Discard supernatant. EXTREMELY CAREFUL!
10. Wash the white pellet with 500 µl 70% EtOH; Centrifuge: 5 min, max, RT
Discard supernatant. EXTREMELY CAREFUL!
11. Vacuum dry: 7-10 min, 35°C, 400-500 mmbar. Make sure no residual EtOH remains. Over-dried pellet may be difficult to resuspend (and difficult to “see”).
12. **Samples on ice!** Completely dissolve the pellet in 30 µl ddH₂O, vortex 2 sec, short spin 2 sec.
13. Nanodrop and check 200 ng on a gel. Over-loading will affect QC assessment.
14. For sequencing purpose, dilute a small aliquot to 200 or 500 ng/µl (depending on plasmid size, refer to sequencing SOP). Store at -20°C (concentrated form is desired for long-term storage).

* Stored at -20°C

* Definitely No PEG! Minor RNA contamination is ok, or go back to #6 (Repeat RNase A)

For 2 samples/cultures:

50 µl	20% SDS
20 µl	10M NaOH
930 µl	s/w